REPRESENTED
山东鹰眼机械科技有限公司
SHANDONG
HAWK-EYE MECHANICAL
TECHNOLOGY CO., LTD.

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Shandong Eagle Eye Machinery Technology Co., Ltd. is located in Yantai, Shandong Province, a national civilized city with beautiful environment and abundant products. The company is mainly engaged in the research and development and production of pharmaceutical machinery and equipment. At present, we have automatic lamp inspection machine series, automatic leak detection machine series, capsule tablet inspection machine, bottle wall cleaning and drying machine, the whole packaging line and other products. Widely used in medicine, food, cosmetics and other industries. The company has obtained: national high-tech enterprises, innovative small and medium-sized enterprises, specialized new small and medium-sized enterprises. Has dozens of invention patent certificates.

The company has advanced production technology, scientific production management system and convenient logistics support. Hawkeye executives are keenly aware of the direction of the market. Joint colleges and universities to develop forward-looking mechanical equipment, so that enterprises more creative, vitality and vitality. Standardization, serialization and generalization are integrated into each of our design drawings; Modern CNC machine tools ensure the precision of parts manufacturing; Close cooperation with major colleges and universities, training high-tech talents, so that Eagle Eye always keep up with the market trend, timely insight into market changes, to achieve customer value creation. Enterprise purpose: "first-class management, first-class technology, first-class products, first-class service" Eagle Eye will be the concept of high quality and low price, good performance, efficient service to relieve you of worries.

The company's main customers are: Sinopsin Group, China Resources Group, Modified Group, Huiren Group, Renhe Group, Kangenbei Group, Broad Group, Harbin Zhenbao Island, Jiangsu Enhua Pharmaceutical, Jimin Trusted Group, Luye Pharmaceutical, Sunflower Pharmaceutical, Shandong Donge Ejiao Group, Tianjin Tisly, Jiangsu Jiuxu Pharmaceutical, Shanghai Baolong Pharmaceutical Guangdong Xiangxue Pharmaceutical, etc.

Welcome new and old customers to visit and guide!

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6S field management Elaborate manufacturing Lean production

机加一角 A corner in mechanical machining

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Hawk Eye Machinery creates safety, stability, and peace of mind for you.

仓库配件 Warehouse accessories

携手高校 Join hands with universities

汇聚全球领先技术 Gather the world's leading technologies

满足客户精准的灯检需求的同时为客户提供整段后道包装线解决方案
While meeting customers' precise lamp inspection needs, we also provide them with complete end-of-line packaging line solutions.

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1-20ml ampoule bottle, Cillin bottle series

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
SD technology: Optical fiber cutting detection technology for visible foreign bodies in pharmaceutical liquid (1/4).

SD can detect movable foreign bodies in liquids with low to medium viscosity and good clarity.

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
The fiber passes through the bottle and is picked up by a sensor at the back.
Movable objects can cut through optical fibers and be captured.
The inspection process is continuous, rather than intermittent photography.

YY-HZDJJ-SD灯检机采用山东鹰眼研发的SD光纤切割光技术检测可动异物,误检率低。
待检产品先被高速旋转,带动起异物旋转。
刹车,瓶子停止转动,而瓶内制品仍旋转。
YY-HZDJJ-SD lamp detector adopts SD light developed by Shandong Eagle Eye
The fiber cutting optical technology can detect movable foreign bodies, and the false detection rate is low.
The product to be inspected is first rotated at high speed to drive the foreign body to rotate.
Brake, the bottle stops rotating, while the contents of the bottle still spin.

检测时间 -与相机对比
Detection time - compared with the camera

连续检测没有间隙
Continuous detection without any gap.
间隙检测存在未检的风险
Gap detection poses the risk of undetected wind. Continuous detection ensures no gaps.

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检测产品类型 APPLICATIONS

1-30ml管制瓶/模制瓶系列
1-30ml tube bottle/molded bottle series

主要性能特点 MAIN CHARACTERS

检测原理 DETECTION PRINCIPLE

机器视觉非接触式检测 Non-contact inspection by machine vision

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(5) 密封性缺陷: 产品密封不良, 瓶内进入空气
- (1) Powder compact defects: Visible foreign substances on the surface, such as glass, lint, black spots, etc.; broken, product discoloration, liquid-like state;

(2) Bottle cap defects: Presence or absence, damage, dirt, scratches, dents of plastic and aluminum caps, missing rubber stoppers, position, etc.;

(3) Bottle body appearance: Cracks, scratches, damage, bubbles, and contaminants inside and outside the container, etc.;

(4) Empty bottles.

(5) Sealing defects: Poor sealing of the product, air entering the bottle.

技术参数 TECHNICAL PARAMETER

外形尺寸(长*宽*高): 1840X2500X2200 (不含进/出瓶网带尺寸, 可按客户要求联线)	External dimensions (length * width * height): 1840X2500X2200 (excluding the dimensions of the bottle inlet and outlet conveyor belts, can be connected as per customer requirements)
净重: ≈3000Kg	Net weight: ≈ 3000 Kg
适用规格: 30ml以内的管制瓶/模制瓶	Applicable specifications: Control bottles/molded bottles of 30 ml or less
检测速度: 600p/min	Detection speed: 600 p/min
电源功率: <12KW	Power supply: < 12 KW
剔除方式: 伯努利气动剔除	Rejecting method: Bernoulli pneumatic rejection
压缩空气: 0.3-0.5MPa	Compressed air: 0.3 - 0.5 MPa
主要材质: 304不锈钢、聚甲醛等	Main materials: 304 stainless steel, polyoxymethylene, etc.
工作温度&湿度: 温度0~35℃, 湿度15%~80%	Working temperature & humidity: Temperature 0 - 35°C, humidity 15% - 80%

规格件为拨轮、栏栅、挡瓶板等, 产品规格变化时需更换规格件。
The specification parts include the dial, the grid, the bottle stopper plate, etc. When the product specifications change, the specification parts need to be replaced.

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This series of lamp detector adopts high-resolution industrial camera, which magnifies the inspected object tens of times, making the visibility of visible foreign objects stronger.

It can detect glass debris, aluminum debris, rubber shaver, hair, fiber and other foreign matter above 40μm in sterile liquid online to ensure the accurate removal of unqualified products, and has independent intellectual property rights of the software patent.

The bottle head is detected by 3 cameras, and the detection effect is better.

Before the bottle enters the testing turntable, the unqualified products of the bottle head are removed to prevent the broken bottles in the testing turntable due to the problems of the bottle head bubble, wire drawing sealing, cracking and so on.

Menu hierarchical system, only the administrator can set parameters to modify, to avoid the operator misoperation to modify the possibility of parameters;

The image of qualified product is saved, which is convenient for users to analyze the existing problems in order to improve the process, and has the function of online detection image playback, which is convenient for online monitoring;

Detection parameters can be freely set, according to customer requirements to achieve the best effect to save; All imported servo motors are used as drivers to ensure the accuracy of equipment operation;

Centrifugal belt adopts well-known brand, with high wear resistance; The operating interface can effectively display the working state of each camera on the device, thus ensuring the safety of product testing;

The electrical control system adopts Omron motion control surface, each motor is independent and accurate output, and the motion position is accurate without misalignment, effectively avoiding the signal delay and delay caused by the traditional PLC triggered signal through the photoelectric sensor, achieving the consistency of the tracking speed and the non-deviation of the millisecond rate of the camera, making the detection effect and equipment stability have a qualitative leap;

Machine fault information feedback setting;

Reserved communication interface for remote fault diagnosis;

Real-time monitoring of light source brightness, camera status, high-speed bottle status, fault information will be displayed on the control panel at the same time alarm shutdown: the height of the machine is completely suitable for general female operators.

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Model:YY-HZDJJLJ(500、600 Model)

Scope of application: 1-20ml ampoules

Test items; foreign matters visible in liquid (such as glass chips, metal chips, hairFiber, white spot, white block), loading capacity, appearance)

Detection speed: 600 bottles / min (taking 2ml as an example)

Resolution:>40 μ m


Micropore resolution:≥0.1μm

Noise≤75db

Power: 9kw

Voltage: AC 380V,50 Hz

Net weight:4200kg

YY-HZDJJ-600X

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AUTOMATIC MULTI-FUNCTIONAL INSPECTION MACHINE

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1. This series of light inspection machines uses high-resolution industrial cameras, magnifying the inspected objects by dozens of times, making visible foreign substances more visible.
2. It can detect foreign substances such as glass shards, aluminum shavings, rubber debris, hair, and fibers larger than 40μm in sterile liquids online, ensuring that substandard products are accurately removed.
3. The software is owned by an enterprise with independent intellectual property rights.
4. The bottle cap and bottle head are inspected by three cameras, providing better inspection results. Unqualified bottle heads are removed before the bottles enter the inspection turntable, preventing broken bottles due to bottle head cracks and affecting the inspection.
5. The menu is hierarchical, and only administrators can set and modify parameters, avoiding the possibility of operators mistakenly modifying parameters.
6. Images of substandard products are saved for users to analyze existing problems and improve the process. It has an online inspection image playback function for convenient online monitoring.
7. Inspection parameters can be freely set and saved when the best effect is achieved according to customer requirements.
8. All imported servo motors are used as the driving force to ensure the accuracy of equipment operation. The centrifugal belt is made of a well-known brand and has extremely high wear resistance.
9. The operation interface can effectively display the working status of each camera on the equipment, ensuring the safety of product inspection.
10. The electrical control system uses Omron motion controllers, with each motor independently and precisely outputting, ensuring accurate movement positions without misalignment. This effectively avoids signal lag and delay caused by traditional PLCs using photoelectric sensors to trigger signals, achieving consistent tracking speed and millisecond-level deviation-free camera image capture, resulting in a significant improvement in inspection results and equipment stability.
11. Machine fault information feedback setting.
12. Communication interfaces are reserved for remote fault diagnosis.
13. Real-time monitoring of light source illuminance, camera status, and centrifugal bottle rotation status. Fault information is immediately displayed on the control panel and an alarm is triggered to stop the machine. The overall height of the machine is completely suitable for female operators.

Mscope of application:1-50 ml
Test items: foreign matters visible in liquid (such as glass chips, metal chips,hair Fiber, white spot, white block), loading capacity,appearance)
Detection speed:400-600 bottles/ min (taking 2ml as an example)
Resolution:>40μm
Noise: 75db
Power:10KW
Voltage:AC 380V,50 Hz
Net weight:4200kg
Owerall dimension:3700x2500x2300 mm(L x W x H]
Table height:930mm

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EYE DROPS AUTOMATIC ROTARY LAMP INSPECTION MACHINE

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1-20ml polyurethane /PVC round bottle eye drops series

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The fully automatic lamp inspection machine employs machine vision technology and automation control principles. It uses mechanical methods to make the inspected products rotate at high speed and brake instantly. It adopts a special light source system and ultra-high-speed precision industrial cameras to obtain a series of images of the inspected products. It uses the intelligent image processing algorithm software system independently developed by the company to compare, calculate and analyze the images, completing the inspection of the products. The substandard products detected are removed by the actuator of the automatic control system.

本系列设备主要由下列部分组成 This series of equipment is mainly composed of the following parts.

- ①进瓶部分:被检产品由进瓶输送带传至进瓶拨盘,再输送至检测旋转台;
- ②检测部分:被检产品高速旋转,经多重特殊光源检测,形成多幅影像供处理系统比对、辨别;
- ③出瓶部分:被检产品经检测后由出瓶拨盘,使用摆臂剔除的方式,将合格品和不合格品分别导入收瓶盘;
- ① Bottle feeding section: The inspected products are conveyed by the bottle feeding conveyor belt to the bottle feeding turntable and then to the inspection rotating platform.
- ② Inspection section: The inspected products rotate at high speed and are inspected by multiple special light sources, forming multiple images for comparison and identification by the processing system.
- ③ Bottle discharging section: After inspection, the inspected products are conveyed by the bottle discharging turntable and the qualified and unqualified products are respectively introduced into the bottle receiving tray by the swing arm rejection method.

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The fully automatic light inspection machine employs machine vision technology and automation control principles. It uses mechanical methods to make the inspected products rotate at high speed and then brake instantly. It adopts a special light source system and ultra-high-speed precision industrial cameras to obtain a series of images of the inspected products. The company's self-developed intelligent image processing algorithm software system is used to compare, calculate and analyze the images to complete the inspection of the products. The unqualified products detected are removed by the actuator of the automatic control system.

Visible foreign matter detection: Based on the principle of machine vision, it goes through a series of processes: "smooth bottle feeding -> high-speed bottle rotation -> emergency stop -> photo taking -> high-speed calculation -> sorting", to separate the products containing visible foreign matter from the normal products.

The interior of the equipment uses a columnar structure, providing a large internal space that is easy to clean and maintain. At the same time, the transparent door design makes it easy to observe the internal operation status.

The slot design makes it convenient and quick to change the specification parts.

The combination of light transmission method and light reflection method and light diffraction method is not smooth. And the special light path setting method makes the detection system of automatic lamp detector stable and efficient.

技术参数 TECHNICAL PARAMETER

规格:3400*2800*2300mm±50mm	Specification: 3400*2800*2300mm ± 50mm
检测项目:玻璃、玻屑、毛发、白点、白块、纤维等可见异物以及口服液瓶盖、包边不良、装量、异色	Inspection items: Visible foreign matters such as glass blocks, glass chips, hair, white spots, white blocks, fibers, etc., as well as oral liquid bottle caps, poor edge wrapping, filling volume, and color differences.
检测速度:600瓶/分钟	Inspection speed: 600 bottles per minute
分辨率:≥40μm	Resolution: ≥ 40μm
噪音:≤75db	Noise: ≤ 75db
功率:9KW	Power: 9KW
光源照度:1000-4000LX	Light source illuminance: 1000 - 4000LX

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检测产品类型 APPLICATIONS

1-50ml西林瓶/口服液瓶系列
1-50ml Cillin bottle/oral liquid bottle series

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
Applicable tablets: plain tablets, film-coated tablets, sugar-coated tablets, single-sided/double-sided printed tablets with score lines and engraved marks

主要性能特点 MAIN CHARACTERS

1. 采用多台CCD多次各象限的重叠检测，漏检率极低。常规要求下检测准确率99.5%。
 2. 可准确的检测出胶囊的各种或内部或外部缺陷，如:叉劈，套环，脏污，小孔，空胶囊，异物胶囊，半填充，顶凹等。
 3. 特别是可以检测印字缺陷。
 4. 运行中漏粉胶囊不会影响到其他胶囊，且非常清洁非常方便。
 5. 搭配真空上料装置可以实现无人操作。
 6. 采用高速高精度的工业CCD及工业电脑，设备运行稳定，检测速度可达200000粒/H，可24小时运转。
1. By using multiple CCDs for repeated quadrant overlap detection, the rate of missed detection is extremely low. Under normal requirements, the detection accuracy rate is 99.5%.
 2. It can accurately detect various internal or external defects of the capsule, such as: split, sleeve ring, dirt, small hole, empty capsule, foreign body capsule, half-filled, top concave, etc.
 3. Especially, it can detect printing defects.
 4. The powder leakage of a capsule during operation will not affect other capsules, and it is very clean and convenient.
 5. When combined with a vacuum feeding device, it can achieve unmanned operation.
 6. Equipped with high-speed and high-precision industrial CCD and industrial computer, the equipment operates stably, with a detection speed of up to 200,000 capsules per hour, and can run for 24 hours.

技术参数 TECHNICAL PARAMETER

设备型号:剂胶囊外观检测机 YY-JPDJJ-01	Model: Capsule appearance testing machine YY-JPDJJ-01
适用剂型:圆形片:直径5-12mm、厚2-8mm异型片:宽5-12mm、厚2-8mm、长5-21mm	Applicable Dosage Forms: Round Tablets: Diameter 5-12mm, Thickness 2-8mm; Irregular Tablets: Width 5-12mm, Thickness 2-8mm, Length 5-21mm
检测项目:污点、刮伤、异物附着、断裂、缺角、变形、异色、变色、印刷模糊、印刷渗漏、包衣不良、刻印不良	Inspection Items: Stains, Scratches, Foreign Matter Adhesion, Breakage, Missing Corners, Deformation, Color Difference, Color Change, Blurred Printing, Printing Leakage, Poor Coating, Poor Engraving
检查面:表背面、侧面	Inspection Surfaces: Front and Back, Sides
检查精度:50μm大小以上的黑点缺陷	Inspection Precision: Black Spot Defects Larger than 50μm
检测速度:5000-10000粒/分钟(以直径5毫米为例)	Inspection Speed: 5000-10000 Pieces/Minute (Take a diameter of 5 mm for example)
灯检亮度:≤ 3000LX。	Light Inspection Brightness: ≤ 3000LX
功率:3KW	Power: 3KW
电压:220V	Voltage: 220V

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6.采用高速高精度的工业CCD及工业电脑，设备运行稳定，检测速度可达200000粒/H，可24小时运转。
1. By using multiple CCDs for repeated quadrant overlap detection, the rate of missed detection is extremely low. Under normal requirements, the detection accuracy rate is 99.5%.

2. It can accurately detect various internal or external defects of the capsule, such as: split, sleeve ring, dirt, small hole, empty capsule, foreign body capsule, half-filled, top concave, etc.

3. Especially, it can detect printing defects.

4. The powder leakage of a capsule during operation will not affect other capsules, and it is very clean and convenient.

5. When combined with a vacuum feeding device, it can achieve unmanned operation.

6. Equipped with high-speed and high-precision industrial CCD and industrial computer, the equipment operates stably, with a detection speed of up to 200,000 capsules per hour, and can run for 24 hours.

技术参数 TECHNICAL PARAMETER

设备型号:剂胶囊外观检测机 YY-JPDJJ-02	Model: Capsule appearance testing machine YY-JPDJJ-02
适用剂型:圆形片:直径5-12mm、厚2-8mm异型片:宽5-12mm、厚2-8mm、长5-21mm	Applicable Dosage Forms: Round Tablets: Diameter 5-12mm, Thickness 2-8mm; Irregular Tablets: Width 5-12mm, Thickness 2-8mm, Length 5-21mm
检测项目:污点、刮伤、异物附着、断裂、缺角、变形、异色、变色、印刷模糊、印刷渗漏、包衣不良、刻印不良	Inspection Items: Stains, Scratches, Foreign Matter Adhesion, Breakage, Missing Corners, Deformation, Color Difference, Color Change, Blurred Printing, Printing Leakage, Poor Coating, Poor Engraving
检查面:表背面、侧面	Inspection Surfaces: Front and Back, Sides
检查精度:50μm大小以上的黑点缺陷	Inspection Precision: Black Spot Defects Larger than 50μm
检测速度:2500-5000粒/分钟(以5毫米为例)	Inspection Speed: 2500-5000 Pieces/Minute (Taking 5mm as an Example)
灯检亮度:≤ 3000LX。	Light Inspection Brightness: ≤ 3000LX
功率:3KW	Power: 3KW
电压:220V	Voltage: 220V

YY-JNJC

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MULTI-CHANNEL CAPSULE INSPECTION MACHINE

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The YY-JNJC7200 capsule appearance inspection machine adopts advanced high-speed cameras for all-round photography. Through AI + CPU processing and algorithmic computation, combined with a stable mechanical structure and a reliable control system, it can detect defects such as top concavity, half capsules, depressions, damages, foreign objects, black spots, hair, abnormal colors, flat shells, fork splits, splits, and wrong insertions of capsules. It conducts high-speed 360-degree precise detection of capsule appearances, quickly identifying defects and accurately and rapidly eliminating unqualified and defective capsules. This solution addresses the issues of subjective emotions, fatigue, accuracy, and speed in manual inspection, significantly reducing the production costs of pharmaceutical enterprises.

设备主要特点 Main features of the equipment:

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⑭ 整个过程无需人工接触胶囊, 减低了人为风险。
- ① The core algorithm is 100% independently developed, with complete independent intellectual property rights.

② Each capsule is inspected 360° without any blind spots, ensuring defective capsules have nowhere to hide.

③ The structure is simple, easy to disassemble and clean; it is convenient to operate and maintain.

④ Compressed air can be utilized to achieve automatic feeding.

⑤ Depending on the number of channels, it can be connected to filling machines of different production capacities or used as a standalone machine.

⑥ It adopts machine vision inspection principles and high-speed image processing technology to automatically screen and remove defective products.

⑦ Inspection information is templated and can be reused; the core inspection software is customized and can be flexibly designed according to the characteristics of the customer's products.

⑧ Inspection data is stored in real time, and inspection results, real-time images, and operational status are displayed dynamically.

⑨ The software is designed in a modular way, making it easy to upgrade, improve, and maintain.

⑩ It can ensure no capsules are left in the channels or feeding hoppers.

⑪ After the capsules are automatically sorted, each one is inspected individually, with a speed of up to 2 capsules per second per channel.

⑫ Applicable capsule sizes: By changing the molds, it can be used for 00#, 0#, 0# extended, 1#, 2#, 3#, 4#, and 5# hard gelatin capsules. Different operation permission levels are set for management and control.

⑬ It is equipped with movable casters for easy and flexible movement.

⑭ The entire process does not require manual contact with the capsules, reducing human risks.

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1. The feeding device, consisting of an in-bottle conveyor belt, a horizontal-to-vertical bottle inverter and an acceleration conveyor belt;
2. The inspection area, where the inspection stations are divided into two parts. The first group of stations adopts an inspection method from bottom to top, mainly inspecting the filling volume, black and white spots on the bottle body and bottom. The second group of stations mainly inspects the appearance of the bottle top, sealing and seams;
3. The bottle ejection section, composed of a reject conveyor belt and an out-bottle conveyor belt. The rejection can achieve classified rejection.

100% comprehensive inspection of BFS ampoules

Reliably identify damage, deformation, foreign particles, incorrect filling, product color deviation, impurities, black burn marks, appearance defects and sealing errors with intelligent LED lighting, optical instruments, image processing and evaluation.

技术参数 TECHNICAL PARAMETER

- ◆稳定产量: 0~150排/分钟 (1ml—10ml)
- ◆Stable output: 0 - 150 rows per minute (1ml - 10ml)
- ◆功率: 4kW, 380V 50Hz
- ◆Power: 4kW, 380V, 50Hz

BFS灯检机设备技术优势 Technical advantages of BFS light inspection machine equipment

- ◆对包装不造成任何损害, 即无损检测。
- ◆Does not cause any damage to the packaging, that is, non-destructive testing.
- ◆对所有产品进行检测, 即在线检测。
- ◆All products are tested, that is, online testing.
- ◆有足够的灵敏度, 能够检测到微小的黑点、白点。(≥40μm异物)
- ◆Has enough sensitivity, can detect small black spots, white spots. (≥40μm foreign body)
- ◆BFS灯检机主要是针对国内联排产品特性, 全方位检测,
- ◆BFS lamp detector is mainly for the characteristics of domestic rooming products, all-round detection,
- ◆检测部位包括: 头部、尾部、正反面。(两组工位, 八个相机)
- ◆Detection parts include: head, tail, front and back. (Two sets of stations, eight cameras)
- ◆能够消除外界因素或人为因素的干扰。
- ◆Can eliminate the interference of external factors or human factors.
- ◆能够经过验证
- ◆Can be verified

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BFS联排塑料安瓿系列
BFS row plastic ampoule series

YY-ADJJ

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AUTOMATIC STRAIGHT LINE LAMP INSPECTION MACHINE

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Freeze-dried powder, powder injection, water injection, oral liquid

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
AUTOMATIC HIGH VOLTAGE DISCHARGE LEAK DETECTION MACHINE

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
The bottle entry and exit adopt the star wheel Bernoulli principle mechanism, with a low bottle breakage rate.
Anti-break bottle device: When a medicine is stuck at the front end, the front end retreats, and the proximity switch senses it, causing the equipment to stop. This prevents bottle breakage.
Safe and reliable rejection system:
Simple parameter setting, users can operate the equipment after a short training period;
Detection electrodes are placed at the bottom, making cleaning convenient. The transformation from upright to lying down and from lying down to upright positions of the medicine is smooth, reducing the risk of bottle breakage;
Equipped with communication interfaces for remote diagnosis;
Detection sensitivity can be adjusted according to customer requirements;
Constantly optimized software system with a user-friendly interface:
Detection sensitivity and parameters can be easily set to meet user needs;
Reserves communication interfaces for remote monitoring and provides free remote online technical services;
Unqualified images are saved, and users can access them online for statistics and process adjustment;
Self-monitoring system during operation (the machine automatically alarms when important operation points fail). The detection speed and sensitivity of the equipment can be adjusted within a certain range as needed;
The equipment uses full servo drive, with stable and reliable actions and basic maintenance-free transmission components;
During the detection process, it can automatically separate defective and qualified products;
Jointly developed high-voltage power supply technology with deep experience and complete intellectual property rights.
Batch detection comparison between human and machine, the equipment's various data are superior to manual detection;
Control system: The equipment uses a PLC system for automatic control, with a 3-level password protection program to prevent parameters from being modified at will.
Mold replacement: The equipment can adapt to different specifications of products by replacing molds. Mold replacement is simple and quick, with accurate positioning, and each specification of molds has clear markings for distinction;
Real-time monitoring of detection voltage. When the voltage is lower than the alarm value, the equipment automatically stops. When the equipment produces consecutive defective products and the quantity exceeds the set value, the equipment automatically stops;
The equipment has a door-opening stop function, and the outer cover is equipped with an electromagnetic protection film;
The equipment has an ozone recovery function;
Emergency stop buttons are installed at the feed port and the bottle collection position at the tail of the operation table;
Convenient for operation and maintenance;
Process data during any alarm will be stored as operation data.

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

检测产品类型 APPLICATIONS

联排塑料安瓿2ml-15ml
Plastic ampoules in a row 2ml - 15ml

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主要性能特点 MAIN CHARACTERS

Shandong Hawk Eye BFS joint leak detection machine is based on Nica, Japan joint plastic ampoule leak detection machine, discharge detection samples through high frequency and high voltage, and distinguish through the current value change of qualified and unqualified products.
Comprehensive detection, mainly detect the head, middle and melt seal. The key technology for row detection accuracy is to detect the corresponding angle and position of the electrodes for each product. Shandong Eagle Eye designs detection electrodes for different shapes of each product. The detection accuracy is more guaranteed.


高压放电检测设备技术优势 Technical advantages of high-voltage discharge testing equipment

- ◆对包装不造成任何损害,即无损检测。
- ◆对所有产品进行检测,即在线检测。
- ◆有足够的灵敏度,能够检测到微小的泄露。(≥3μm100%检出)
- ◆BFS检漏机主要是针对国内联排产品特性,全方位检测,检测部位包括:头部、尾部、正反面。(4电极全覆盖)
- ◆能够消除外界因素或人为因素的干扰。
- ◆能够经过验证。
- ◆No damage is caused to the packaging, namely, nondestructive testing.
- ◆All products are tested, that is, online testing.
- ◆With sufficient sensitivity to detect minor leaks.(3 μ m100% detected)
- ◆BFS leak detection machine is mainly aimed at the characteristics of domestic row products, all-round detection, detection parts include: head, tail, front and back.(4 electrode full coverage)
- ◆Eliminating the interference of external factors or human factors.
- ◆Be verified.

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- ◆Both positive and negative electrodes shall contact drugs or have very close distance (2-3mm) and the distance is adjustable; to facilitate stable detection loop during use.
- ◆There is enough distance between the positive and negative poles to ensure that the positive and negative poles are not directly connected when testing.

技术参数 TECHNICAL PARAMETER

设备型号:BFS 高压放电检漏机	Backup model :BFS high voltage discharge leak detector
稳定产量:0~140排/分钟	Stable output: 0~140 rows/min
电容量: 2kW (含高频高压电源功率), 380V 50Hz	Capacity: 2kW (including high frequency and high voltage power supply), 380V 50Hz
高频:600Hz	High frequency :600Hz
高压:0~22kV可调	High voltage :0~22kV adjustable
高频高压电源功率:200W	High frequency high voltage power supply :200W

YY-PTJCJ

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AUTOMATIC BOTTLE HEAD TESTING MACHINE

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
1-20ml ampoules, vials / oral liquid bottles series

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The cleaning and drying process is divided into four stations:

Rough washing station: the bottle enters the bottle from the bottle inlet mesh belt device into the bottle into the bottle large dial, the bottle in the process of progress, to the designed washing station, 4 directions of warm water (temperature adjustable) pre-rough washing, in the rough washing process, the installation of a brush for all-round cleaning of the bottle outside the wall.

Fine washing station: The bottle from the bottle tray into the transition bottle tray, the bottle in the process of progress, to the design of the washing station, 4 directions purified water fine washing. Water washing station adopts pure water high pressure cleaning, and uses the high pressure of the pump to wash the outer wall and the bottom of the bottle. The spray nozzle has a unique shape design, which can wash the outer wall and the bottom of the bottle in all directions. Rinse the water falling from the bottle, recycle it to the rough washing tank, reduce energy consumption, the water in the rough washing tank, through the overflow pipeline, discharge to the floor drain.

Hot air blowing station: The bottle comes out of the water fine washing station and enters the hot air blowing station. The hot air blowing station adopts the company's unique air heating technology, and uses the hot air to blow the water attached to the outside wall of the bottle. The air nozzle adopts a unique shape design, which can blow the water on the outer wall of the bottle in an all-round and multi-angle.

Compressed air blowing station: the bottle from the transition bottle tray is transferred to the bottle out of the bottle pan, into the compressed air blowing station, compressed air uses clean compressed air to blow the water attached to the outside wall of the bottle and the bottom of the bottle. The air nozzle has a unique shape design, which can blow the water from the outer wall and the bottom of the bottle in an all-round way.

Bottle out process: The bottle is completed by air blowing and sent to the cleaning machine through the bottle out track. Can be connected to the next process.

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Equipment model :YY-PWQX
Stable production speed :600 bottles/min
Total power :7.2Kw 380V/50HZ380V/50HZ three-phase five-wire (Depending on customer needs, water heating can be turned off)
Compressed air :0.6Mpa 700L/min
Pure water :0.2~ 0.3Mpa 1000L/h
Weight :1500Kg

Represented by Harikrushna Machines Pvt. Ltd.

+91 98795 44851 marketing@harikrushna.com www.hmplsterile.com | www.harikrushna.com

Plot No. 513, Phase IV, Vatva GIDC, Ahmedabad, Gujarat 382445

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SHANDONG HAWK-EYE MECHANICAL TECHNOLOGY CO.,LTD.

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Add: No. 157, Longchang Road, Laishan District, Yantai City, Shandong Province

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网址website: www.sdyydjj.com